



## R.O. SYSTEM SERVICE LOG

### FACILITY:

|                       |                           |               |              |  |
|-----------------------|---------------------------|---------------|--------------|--|
| <b>SERIAL #</b>       |                           | <b>BEFORE</b> | <b>AFTER</b> |  |
| <b>MODEL #</b>        |                           |               |              |  |
| <b>DATE</b>           |                           |               |              |  |
| <b>SERVICE REPORT</b> |                           |               |              |  |
|                       |                           |               |              |  |
| PRE-FILTER            | (MIN. 20 PSI, MAX 80 PSI) |               |              |  |
| POST FILTER           | (MIN. 20 PSI, MAX 80 PSI) |               |              |  |
| FEED TDS.             | CH. 1                     |               |              |  |
| PROD. TDS.            | (MAX. 100 RPM) CH. 2      |               |              |  |
| % REJECT              |                           |               |              |  |
| MOD. IN PRES.         | (MAX. 350 PSI)            |               |              |  |
| MOD. OUT PRES.        | (MAX. 350 PSI)            |               |              |  |
| REJECT FLOW           | GPM MIN.                  |               |              |  |
| PRODUCT FLOW          | GPM MIN.                  |               |              |  |
| WATER TEMPERATURE     | (MAX 38 C)                |               |              |  |
| STERILIZED            |                           |               |              |  |
| PYROGEN RESULTS       |                           |               |              |  |
| BACTERIA RESULTS      |                           |               |              |  |
|                       |                           |               | TECH. INIT.  |  |

LUBRICATE TONKAFLO PUMP EVERY 4400 HOURS