



R.O.SYSTEM DAILY LOG

FACILITY:

YEAR:		MON.	TUES.	WED.	THURS.	FRI.	SAT.	SUN.
DATE:								
HOUR METER READING								
A.M. HARDNESS TEST (1 GR. MAX.)								
P.M. HARDNESS TEST (1 GR. MAX.)								
SHIFT 1 (MAX. 0.1 PPM) CHL.								
TIME								
INIT.								
SHIFT 2 (MAX. 0.1 PPM) CHL.								
TIME								
INIT.								
SHIFT 3 (MAX. 0.1 PPM) CHL.								
TIME								
INIT.								
SHIFT 4 (MAX. 0.1 PPM) CHL.								
TIME								
INIT.								
PRE-FILTER (MIN. 20 PSI, MAX 80 PSI)								
POST FILTER (MIN. 20 PSI, MAX 80 PSI)								
FEED TDS. CH. 1								
PROD. TDS. (MAX. 100 PPM) CH. 2								
% REJECT								
MOD. IN PRES. (MAX. 350 PSI)								
MOD. OUT PRES. (MAX. 350 PSI)								
REJECT FLOW GPM MIN.								
PRODUCT FLOW GPM MIN.								
PUMP LUBE								
WATER TEMPERATURE (MAX 38 C)								
STERILIZED								
PYROGEN RESULTS								
BACTERIA RESULTS								
TECH. INIT.								

data/rosdl

LUBRICATE TONKAFLO PUMP EVERY 4400 HOURS



(800) 755-1707